

CARVENE

GUTTA PERCHA SOLVENT

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Product Description :

Carvene Gutta Percha Solvent is used for removing gutta percha at the time of root canal re-treatment. Also useful for removing all types of root canal filling materials including Zinc Oxide – Eugenol Cement or as a cleaner for endodontic instruments. Carvene is non-hazardous unlike Chloroform/Trichloroethylene which are considered as carcinogenic substances. Also it is less volatile as compared to chloroform / Trichloroethylene. Carvene Gutta Percha Solvent is based on natural substances and therefore it does not cause irritation to oral mucosa and peri-apical tissue.

Composition : d-Limonene which is a natural substance extracted from orange oil.

Indications :

Carvene G.P. Solvent is used for removing gutta percha at the time of root canal re-treatment. Also useful for removing all types of root canal filling materials including Zinc Oxide – Eugenol Cement or as a cleaner for endodontic instruments.

Instructions for Use : Dispense Carvene Gutta Percha Solvent into dappen dish or directly on cotton pallet. Place the soaked cotton pellet on root canal filling material and cover with temporary cement (Zinc Phosphate). Remove the softened filling material after two days.

Warning : Avoid prolonged or repeated contact with skin (allergic contact dermatitis), and eyes. Avoid prolonged inhalation. Do not take internally Contact with the skin may cause irritation, burns of hypersensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

Cautions : Carvene Gutta Percha Solvent is volatile and flammable. Keep away from heat source. Keep bottle tightly capped after every use.

Storage Conditions : Store at temperature between 10°C to 24°C. Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place.

Shelf life : 3 years from the date of manufacturing.

Presentation :

Ref : 40030

1 x 10ml Glass Bottle, 1 x Dispensing Dropper.

Disposal conditions : To be disposed off based on the prevailing rules of the country in which it is used.

 **Prevest DenPro Limited**
Unit II, EPIP, Bari Brahmna,
Jammu-181133, India.
www.prevestdenpro.com
Customer Care No. : 1800-180-7172
Email : customercare@prevestdenpro.com




Obelis s.a
Boulevard General Wahnis 53
1030 Brussels, Belgium
Tel: +(32) 2 732 99 54, Fax: +(32) 2 732 60 03
E-mail : mail@obelis.net


0123

Revision :	05
Date :	23-05-2019

 Symbol for "BATCH CODE"  Symbol for "USE BY"  Symbol for "NON STERILE"  Symbol for "KEEP DRY"

 Symbol for "CAUTION"  Symbol for "CONSULT INSTRUCTION FOR USE"  Symbol for "CATALOGUE NUMBER"

 SYMBOL FOR EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE  Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"

 Symbol for "KEEP OUT OF THE REACH OF CHILDREN"  Symbol for "MANUFACTURER NAME & ADDRESS"

 Symbol for "TEMPERATURE LIMITATION"  Symbol for "KEEP AWAY FROM SUNLIGHT"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.