

EpoSolv

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Solvent For Resin & Phenolic Compounds

Product Description : EpoSolv solvent is used for removing for resin sealers and phenolic compounds at the time of root canal re-treatment.

Composition : Formamide and Phenylethanol.

Indications : EpoSolv is a solvent for epoxy and other resin sealers during endodontic re-treatment. EpoSolv aids in the softening of phenolic resin type fillings. EpoSolv does not provoke excessive salivation if it comes into contact with the mucosa.

Instructions for Use :

Recommended one visit method :

Clean the pulp-chamber and entrance to the canal. Place a drop of EpoSolv into the chamber using a micro applicator tip. Entry into the canal should be done by once dipping the point of the instrument in EpoSolv before each application. The instrument must be removed from the canal as soon as a certain resistance is felt, then wipe and repeat the operation a number of times, until reaching the apex.

Recommended two visit method : (For the removal of the root canal filling material at a subsequent visit).

On the initial visit : By mechanical means, clean the pulp-chamber and entrance to the canal. Place into the cavity a small cotton pellet soaked in EpoSolv. Compress slightly. Fill with a sealing cement. Re-open either the next day or some days later. The paste will have disintegrated and the mechanical reaming will present no difficulty.

Warning : Avoid prolonged or repeated contact with skin (allergic contact dermatitis), and eyes. Avoid prolonged inhalation. Do not take internally. Contact with the skin may cause irritation, burns of hyper sensitivity. If contact with skin occurs immediately wipe off thoroughly with cotton and alcohol and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention. If the product comes into contact with the eyes, rinse thoroughly with fresh water and consult a physician. It should not be used with patients who have a history of severe allergic reaction to any of the components.

Cautions :

A tooth treated with EpoSolv must be carefully dried before a new restoration made of eugenol resin or acrylic, is placed.

Storage Conditions :

Store at temperature between 10°C to 24°C. Close the cap carefully after use. Keep away from moisture. Shelf life not affected after opening if stored in a cool dry place.

Shelf life : 3 years from the date of manufacturing.

Presentation :

Ref : 40035

1 x 15ml Glass Bottle, 1 x Dispensing Dropper.

Disposal Conditions : To be disposed off based on the prevailing rules of the country in which it is used.















Prevest DenPro Limited
Unit II, EPIP, Bari Brahmna,
Jammu-181133, India.
www.prevestdenpro.com
Customer Care No. : 1800-180-7172
Email : customercare@prevestdenpro.com

EC REP
Obelis s.a
Boulevard General Wahis 53
1030 Brussels, Belgium
Tel: +(32) 2 732 59 54, Fax: +(32) 2 732 60 03
E-mail : mail@obelis.net

 Symbol for "HARMFUL IRRITANT"


0123

Revision :	04
Date :	23-05-2019

-  Symbol for "BATCH CODE"  Symbol for "USE BY"  Symbol for "NON STERILE"  Symbol for "KEEP DRY"
-  Symbol for "CAUTION"  Symbol for "CONSULT INSTRUCTION FOR USE"  Symbol for "CATALOGUE NUMBER"
-  SYMBOL FOR EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE  Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"
-  Symbol for "KEEP OUT OF THE REACH OF CHILDREN"  Symbol for "MANUFACTURER NAME & ADDRESS"
-  Symbol for "TEMPERATURE LIMITATION"  Symbol for "KEEP AWAY FROM SUNLIGHT"  Symbol for "DANGER"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.