

# ViroGuard

**PREVESTDenPro®**  
THE FUTURE OF DENTISTRY

## Surface Disinfectant Spray for Medical Devices

**Product Description :** ViroGuard is a ready to use fast acting and broad spectrum spray for quick cleaning and disinfecting of medical devices and surfaces after each session in medical and dental offices & clinics. It is specially designed to avoid cross contamination in between patients.

### Composition :

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2-Butoxyethanol Surfactant, IPA, pH modifier and water.

**Intended Purpose/Intended Use:** It is intended to be used as surface spray for cleaning, decontaminating and disinfecting medical device and surfaces.

### Indications for use:

It is used for cleaning, decontaminating and disinfecting surfaces (such as trolleys, beds, chairs, trays, work benches, examination tables, dental units, handles, spittoons, other non invasive equipments and instruments) in direct contact with patients and medical personnel.

**Instructions for Use :** Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 sec for effective disinfection and wipe with the clean lint free disposable cloth. (Note: Do not rinse).

**Warnings :** Harmful if swallowed. Causes severe skin burns and eye damage.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF in EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and continue rinsing. Get medical advice/attention if you feel unwell

### Precautions :

- ♦ Advised for use on surfaces of medical devices by professionals only.
- ♦ Wear personal protective equipment such as gloves, masks, goggles for protection. Wash your hands after handling.
- ♦ Do not mix with any other disinfectants.

**Contraindications :** Avoid heat

**Storage Conditions :** Store at temperature between 10°C(50°F) to 24°C(75°F). Keep away from moisture and heat. Close the cap carefully after use. The shelf life of the product is valid only if left unopened and correctly stored in a cool dry place.

**Shelf life :** 3 years from the date of manufacturing.

### Presentation :

Ref: 16012

500ml Spray Bottle

Ref: 16013

1Litre Can

Ref: 16014

5 Litre Can

**Disposal Conditions :** To be disposed off based on the prevailing rules of the country in which it is used.

### Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,  
(EPIP), Bari Brahmna, Jammu-181133,  
India. [www.prevestdenpro.com](http://www.prevestdenpro.com)








Customer Care No. : + 91 8899074151

Email : [customercare@prevestdenpro.com](mailto:customercare@prevestdenpro.com)  
Made in India

CE

Revision :	00
Date :	06-05-2023

**UDI** Symbol for "UNIQUE DEVICE IDENTIFIER"

- LOT** Symbol for "BATCH CODE"  Symbol for "USE-BY DATE"  Symbol for "NON STERILE"  Symbol for "KEEP DRY"
-  Symbol for "CAUTION"  Symbol for "CONSULT INSTRUCTION FOR USE" **REF** Symbol for "CATALOGUE NUMBER"
- EC REP** Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"
-  Symbol for "DO NOT USE IF PACKAGE IS DAMAGED"  Symbol for "KEEP OUT OF THE REACH OF CHILDREN"  Symbol for "MANUFACTURER"
-  Symbol for "TEMPERATURE LIMIT"  Symbol for "KEEP AWAY FROM SUNLIGHT" **MD** Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.