

ViroGuard Foam

PREVESTDenPro®
THE FUTURE OF DENTISTRY

Alcohol Free Disinfectant Foam Spray For Medical Devices

Product Description : ViroGuard Foam is an alcohol free formation for rapid disinfection of alcohol sensitive non invasive medical devices and medical equipments having surface made of glass, rubber, acrylic plastic & artificial leather upholstery.

Composition :

N-(3-aminopropyl)-N-dodecylpropano-1, 3-Diamine, Didecyl-dimethyl Ammonium Chloride, 2-Butoxyethanol Surfactant, D-Glucopyranose, oligomers, decyl octyl glycosides, Sodium lauroyl sarcosinate, pH modifier and water

Intended Purpose/Intended Use: It is intended to be used as a surface spray for cleaning, decontaminating and disinfecting medical device and surfaces.

Indications for use:

- ◆ It is an alcohol free foam spray used for disinfection of alcohol sensitive non invasive medical devices such as dental units, beds, stretches, work benches etc.
- ◆ It can also be used for disinfection of non invasive medical equipments such as monitors, displays, keyboards and control panels etc.

Instructions for Use : Spray directly on the surfaces and objects to be disinfected and leave the surfaces wet for 30 sec for effective disinfection and wipe with the clean lint free disposable cloth. For difficult surfaces such as keyboards, displays, monitors and control panels, use wet clean disposable cloth with spray foam and thoroughly wipe the surfaces.

Warnings : Harmful if swallowed. Causes severe skin burns and eye damage. May cause damage to organs through prolonged or repeated exposure.

IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

IF ON SKIN(or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and continue rinsing. Get medical advice/attention if you feel unwell

Precautions :

- ◆ Advised for use on surfaces of medical devices by professionals only.
- ◆ Wear protective gloves and garments and protect the eyes and face. Wash your hands after handling.
- ◆ Do not rinse or mix with any other disinfectants.

Contraindications : Avoid heat

Storage Conditions : Store at temperature between 10°C(50°F) to 24°C(75°F). Keep away from moisture and heat. Close the cap carefully after use. The shelf life of the product is valid only if left unopened and correctly stored in a cool dry place.

Shelf life : 3 years from the date of manufacturing.

Presentation :

Ref : 16009

500ml Spray Bottle

Ref : 16010

1 Litre Can

Ref : 16011

5 Litre Can

Disposal Conditions : To be disposed off based on the prevailing rules of the country in which it is used.

Prevest DenPro Limited

Unit II, Export Promotion Industrial Park,
(EPIP), Bari Brahmna, Jammu-181133,
India. www.prevestdenpro.com

Customer Care No. : + 91 8899074151

Email : customercare@prevestdenpro.com
Made in India

CE

Revision :	00
Date :	25-04-2023

 UDI Symbol for "UNIQUE DEVICE IDENTIFIER"

 LOT Symbol for "BATCH CODE"  Symbol for "USE-BY DATE"  Symbol for "NON STERILE"  Symbol for "KEEP DRY"

 Symbol for "CAUTION"  Symbol for "CONSULT INSTRUCTION FOR USE"  REF Symbol for "CATALOGUE NUMBER"

 EC/REP Symbol for "AUTHORIZED REPRESENTATIVE IN THE EUROPEAN COMMUNITY/EUROPEAN UNION"

 Symbol for "DO NOT USE IF PACKAGE IS DAMAGED"  Symbol for "KEEP OUT OF THE REACH OF CHILDREN"  Symbol for "MANUFACTURER"

 Symbol for "TEMPERATURE LIMIT"  Symbol for "KEEP AWAY FROM SUNLIGHT"  MD Symbol for "MEDICAL DEVICE"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding. For feedback, complaints and reporting of unexpected events contact customer care.